

Aim: To determine the amount of Dissolved oxygen in given water sample by winkler's method

Apparatus Required: Burette, Pipette, Conical flask, Cork, Beaker, Volumetric flask, Glass rod, Weighing bottle

Chemicals Required:

1. Manganous Sulphate Solution ($\text{MnSO}_4 \cdot 4\text{H}_2\text{O}$)
2. Alkali-iodine Oxide
3. Conc. Sulphuric Acid (H_2SO_4)
4. Starch Solution
5. Sodium Thiosulphate Solution (0.025 N)
6. Potassium Dichromate Solution (0.025 N)

Dissolved Oxygen

- ❖ Oxygen is one of the most common dissolved gases in the water
- ❖ DO is the most important indicator of the health of a water body
- ❖ DO is vital for survival of aquatic life in water bodies.
- ❖ A higher dissolved oxygen level indicates a better water quality. If dissolved oxygen levels are too low, some fish and other organisms may not be able to survive.
- ❖ Oxygen affects a vast number of other water indicators, not only biochemical but also like the odor, clarity and taste.

Dissolved Oxygen

- ❖ D.O. range 5-9 mg/L: Good quality water
- ❖ Value below 5mg/L: Aquatic organism become stressed
- ❖ Value below 2mg/L: hypoxic water cause "Mass Fish Kill"


Factor Affecting Dissolved Oxygen

- ❖ Temperature
- ❖ Pressure
- ❖ Salinity
- ❖ Surface area
- ❖ Level of organic activities

Oxygen can become dissolved in three ways

- ❖ **Introduced into water by algae, through photosynthesis**
- ❖ **Enters water directly from atmosphere**
- ❖ **Introduced by mechanical equipment**

Dissolved Oxygen Measurement Methods



Titration method/ Winkler method/ Iodometric method for determination of Dissolved Oxygen

- It is not possible to directly measure the amount of dissolved oxygen in a water sample directly.
- The dissolved oxygen does not directly react with another suitable reagent, an indirect procedure was developed by Winkler.
- An iodine/thiosulfate titration can be used to measure the dissolved oxygen present in a water sample.

Titration method/ Winkler method/ Iodometric method for determination of Dissolved Oxygen

- 1. Preparation of standard solution of $K_2Cr_2O_7$ sample**
- 2. Standardization of $Na_2S_2O_3$**
- 3. Preparation of water sample**
- 4. Titration of water sample against $Na_2S_2O_3$**

1. Preparation of standard solution of $K_2Cr_2O_7$ sample (0.02 N)

2. Standardization of $Na_2S_2O_3$:

Take 10 mL standard $K_2Cr_2O_7$ solution in the conical flask + Conc H_2SO_4 (5 mL) + 5% KI (5 mL).

↓ Keep in dark for 15 minutes

Brown Ppt

↓
Titrate with $Na_2S_2O_3$

↓
Pale Yellow Colour

↓ Add 2-3 drops of starch,

↓ Titrate with $Na_2S_2O_3$

Colorless

↓
Repeat it till 3-concordant reading

3. Preparation of water sample :

Fill the reagent bottle with water

↓
Remove 5 mL water by pipette+1 mL
 MnSO_4 + 1 mL Alkali Iodide azide

Replace Stopper & invert several times to mix

↓
Flock in the solution is
allowed to settle

After settling add 1 mL Conc. H_2SO_4

↓
Stopped it & invert several times to mix
& dissolved all the flock (Repeat it till
ppt dissolve)

4. Titration of water sample against $\text{Na}_2\text{S}_2\text{O}_3$:

Take 100 mL water sample



Titrate with $\text{Na}_2\text{S}_2\text{O}_3$



Yellow Colour



Add 2- drops of starch,
Titrate with $\text{Na}_2\text{S}_2\text{O}_3$


Colorless solution




Repeat it till 3-concordant
reading

Chemical Reactions

Manganese sulfate in alkaline conditions



Mix with sample under water – This reacts with the dissolved oxygen to produce a brown precipitate.



Adding concentrated H_2SO_4 – enables the Mn(IV) compound to release free iodine from KI.



Titration : The free iodine is then titrated with standard sodium thiosulfate




Table -1 Standardization of $\text{Na}_2\text{S}_2\text{O}_3$ using starch as an indicator

S. No.	Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ (mL)	Burette reading (mL)		Titrant used Difference (mL)
		Initial	Final	
1.	10			
2.	10			
3.	10			
4.	10			
5.	10			

Table -2 Titration of water sample with $\text{Na}_2\text{S}_2\text{O}_3$ using starch as an indicator

S. No.	Volume of Water sample (mL)	Burette reading (mL)		Titrant used Difference (mL)
		Initial	Final	
1.	100			
2.	100			

$\text{MnSO}_4 \rightarrow$ Catalyst for the R_2O_7
alkali azide \rightarrow

$$N_1 (\text{K}_2\text{Cr}_2\text{O}_7) = \frac{76}{\text{eqwt}} \times \frac{1000}{100}$$

=

$$N_1 V_1 = N_2 V_2$$

$(\text{K}_2\text{Cr}_2\text{O}_7)$ $(\text{Na}_2\text{S}_2\text{O}_3)$

$$N_1 \times 10 = N_2 \times \text{Vol}^{\text{m}} \text{ consume}$$

$$\begin{aligned} &= N_2 \\ &= N_2 V_2 \\ &= \text{water sample} \end{aligned}$$

$$N_2 \times \text{Vol}^{\text{m}} \text{ consume} = N_3 \times 100 \text{ ml}$$

$$= N_3$$

Amount of dissolved oxygen = $N_3 \times 1000 \times 8$

$$= N_2 \times 1000 \times 8 = \text{mg/l}$$

N_1 = Normality of $\text{K}_2\text{Cr}_2\text{O}_7$

V_1 = Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ in conical flask

N_2 = Normality of $\text{Na}_2\text{S}_2\text{O}_3$

V_2 = Volume of $\text{Na}_2\text{S}_2\text{O}_3$ consumed (from first observation table)

N_3 = Normality of water sample

V_3 = Volume of water sample

V_2 = Volume of $\text{Na}_2\text{S}_2\text{O}_3$ consumed (from second observation table)

RESULT

The Dissolved oxygen of given water sample is found to be ----- mg/L