






Aim: To estimate the amount of Fe(II) ions present in the unknown solution of Mohr's salt by titrating it with $\text{K}_2\text{Cr}_2\text{O}_7$ ($\sim \text{N}/40$) using internal indicator.

Requirements:

Apparatus: Burette, pipette, conical flask, standard flask, beaker, spatula, wash bottle, weighing bottle.

Chemicals: Mohr's salt solution, $\text{K}_2\text{Cr}_2\text{O}_7$ solution ($\sim \text{N}/40$), o-phosphoric acid, dilute sulphuric acid, Diphenyl amine indicator and Mohr's salt (solid).



Preparation of $\left(\frac{N}{40}\right)$ Mohr's Salt Solution (100ml)



$$\text{Normality} = \frac{\text{wt.} \times 1000}{\text{Eq. wt.} \times \text{V(ml)}}$$

$$\frac{1}{40} = \frac{\text{wt.} \times 1000}{392.14 \times 100}$$

$$\text{wt.} = \frac{392.14}{400}$$

$$\boxed{\text{wt.} = 0.98 \text{ g}}$$

Equivalent wt. of $\text{K}_2\text{Cr}_2\text{O}_7$:





$$\text{Molecular wt. of } \text{K}_2\text{Cr}_2\text{O}_7 = 294.21$$

$$\begin{aligned} \text{Equivalent wt. of } \text{K}_2\text{Cr}_2\text{O}_7 &= \frac{\text{Molecular wt.}}{\text{change in o.s.}} \\ &= \frac{294.21}{6} \end{aligned}$$

$$\boxed{\text{Equivalent wt. of } \text{K}_2\text{Cr}_2\text{O}_7 = 49.03}$$

Role of Phosphoric Acid:



2. H_3PO_4 also prevents the premature oxidation of the indicator

Procedure:

1. Preparation of Standard Solution of Mohr's salt:

Prepare a standard solution of Mohr's salt $\frac{N}{40}$ by dissolving 0.98g of Mohr's salt in $\sim 10\text{ mL}$ of dil. H_2SO_4 then dilute with water to make 100 mL of solution.

2. Standardisation of Potassium Dichromate:

Pipette out 10 mL of standard Mohr's salt in a conical flask. Add 10 mL dil H_2SO_4 in the Mohr's salt solution then add 2 mL H_3PO_4 and diphenyl amine (2-3 drops), respectively. Titrate with $\text{K}_2\text{Cr}_2\text{O}_7$ solution till a permanent blue colour is obtained. Repeat the titration to get three concordant values.

Observation Table:

Sr.No.	Volume of Mohr's salt (mL)	Burette Reading (mL)		Volume of $\text{K}_2\text{Cr}_2\text{O}_7$ used (mL)
		Initial	Final	
1.	10.0	0.0	9.1	10.3 mL
2.	10.0	9.1	19.4	
3.	10.0	19.4	29.7	
4.	10.0	29.7	40.0	

Concordant volume = 10.3 mL

$$N_1 V_1 = N_2 V_2$$

$$\text{K}_2\text{Cr}_2\text{O}_7 = \text{Standard Mohr's salt solution}$$

$$N_1 \times 10.3 = \frac{1}{40} \times 10$$

$$N_1 = \frac{1}{41.2} N$$

$$\text{Normality of } \text{K}_2\text{Cr}_2\text{O}_7 \text{ solution} = \frac{1}{41.2} N$$

3. Titration of Unknown Mohr's salt Solution with Standard $K_2Cr_2O_7$ Solution:

Pipette out the unknown Mohr's salt solution (10 mL) in a conical flask. Add dil. H_2SO_4 (10 mL) & then add 2 mL H_3PO_4 and 2-3 drops diphenyl amine (2-3 drops) indicator. Titrate with standard $K_2Cr_2O_7$ solution till a permanent blue colour is obtained. Repeat the titration to get three concordant values.

Observation Table 2:

Sr.No.	Volume of Mohr's Salt (mL)	Burette Reading (mL)		Volume of $K_2Cr_2O_7$ Used (mL)
		Initial	Final	
1.	10.0	0.0	11.5	12.0 mL
2.	10.0	11.5	23.5	
3.	10.0	23.5	35.5	
4.	10.0	35.5	47.5	

Concordant Volume = 12 mL

$$\begin{aligned}
 N_1 V_1 &= N_2 V_2 \\
 \text{Standard } K_2Cr_2O_7 &= \text{Unknown Mohr's Salt Solution} \\
 \frac{1}{41.3} \times 12 &= N_2 \times 10
 \end{aligned}$$

$$N_2 = \frac{12}{413}$$

$$N_2 = 0.029 \text{ N}$$

Normality of unknown Mohr's salt solution = 0.029 N
Calculations:


$$\text{Strength of Mohr's salt solution} = \text{Normality} \times \text{Equivalent wt. of Mohr's salt}$$

$$\text{Strength in g/L} = 0.029 \times 392.14$$

$$= 11.393 \text{ g/L}$$

$$\text{Amount of Mohr's salt present} = 11.393 \text{ g/L}$$

Amount of Fe(II):



$$\text{Amount of Fe(II) in 1 g Mohr's salt} = \frac{56}{392.14}$$

$\therefore$ Amount of Fe(II) present in 11.393 g Mohr's salt

$$= \frac{56}{392.14} \times 11.393$$

$$\text{Amount of Fe(II)} = 1.626 \text{ g}$$

Result: Amount of Fe(II) present = 1.626 g