

Date _____
Page _____




Aim :- To determine the strength of unknown oxalic acid by titrating with Potassium Permanganate solution. Prepare standard oxalic acid solution of $\frac{N}{40}$

Apparatus Required: Burette, Pipette, weighing bottle, washing bottle, iron stand, Beaker - 400ml. (2) wire gauze, tripod stand, conical flask.

Theory $\rightarrow$ It is reduction-oxidation titration where Potassium Permanganate gain $5e^-$ while oxalic acid undergoes oxidation give 2 electrons. $KMnO_4$ undergoes reduction.



The reducing agent is -----

The oxidizing agent is -----






Since nascent (O) is not balanced in eq (1) & (2) therefore we multiply eq (2) by (5)

Ionic Equation




Multiply eq. (1) by 2 and eq. (2) by 5




Procedure

Preparation of standard oxalic acid
standard flask, weighing bottle, spatula funnel
wash bottle is required




weighing bottle



standard flask



Funnel



weighing balance

Since $\frac{N}{40}$ oxalic acid standard solution is to be prepared.

Therefore calculation of mass of oxalic acid is required for its preparation.

Molecular weight of oxalic Acid = 126 g/mole

To prepare $\frac{N}{40}$ normal solution, equivalent weight of oxalic acid is required.

In the eq. (2) of O.A. there is $2e^-$

Date _____
Page _____

are released hence to calculate its equivalent weight its molecular weight is divided by that number

$$\begin{aligned}\text{Equivalent weight} &= \frac{\text{Molecular weight}}{2} \\ &= \frac{126 \text{ g/mole}}{2} \\ &= 63\end{aligned}$$

To prepare 1 Normal solution i.e. 63 g of oxalic acid is dissolved in 1000 ml of H_2O (1 l H_2O)
For 1000 ml.

$$0.1 \text{ N} \longrightarrow 6.3 \text{ g}$$

$$0.01 \text{ N} \longrightarrow 0.63 \text{ g}$$

$$\frac{\text{N}}{40} \longrightarrow \frac{63}{40} = 1.575$$

For $\frac{\text{N}}{40}$ oxalic acid in 1000 ml $\rightarrow 1.575 \text{ g}$ is to be dissolved 1000 ml of H_2O
Since we have to do maximum 3 to 5 titrations with this solution hence there is no requirement of 1000 ml solution. 100 ml of standard solution is sufficient


$$1000 \text{ ml} \rightarrow 1.575 \text{ g}$$

$$100 \text{ ml} \rightarrow \frac{1.575}{1000} \times 100 = 0.1575 \text{ g}$$

First mass of empty weighing bottle is done. Let us suppose it is 5g = m,
Then in the weighing bottle 0.1575g of oxalic is added with help of spatula & its weighed is done 5.1575g

Date _____
Page _____

This mass is transferred to the standard flask with help of funnel.



After transferring the empty weighing bottle is again weighed to know the mass of weighing bottle after transference so that exact amount of oxalic acid transferred can be known.


Let the mass of weighing bottle after transference is 5.0015 g .

$$\therefore \text{Mass of oxalic acid transferred in standard flask is} = 5.1575 \text{ g} - 5.0015 \text{ g} \\ = 0.1560 \text{ g}$$

This amt 0.1560 g of oxalic acid is transferred to standard flask.

Add 10 ml of H_2O from the funnel so that some particles of oxalic acid which are attached in the neck of funnel can be removed & go into standard flask.

Dissolve O.A. in minimum amt of H_2O and then make the solution upto mark.



Normality of the standard solution is



$$\frac{\text{Wt. of O.A.} \times 1000}{\text{Equivalent wt} \times 100}$$

$$= \frac{0.156 \times 1000}{63 \times 100}$$


$$= \frac{1.56}{63} = 0.02476 \text{ N}$$

Step 2. Titration of Potassium Permanganate VS Standard Oxalic Acid solution

In the burette Potassium Permanganate (Prepared by lab staff) is added.



Pipette out 10 ml of standard oxalic acid solution whose normality is 0.02476 and add half T.T. of H_2SO_4 sulphuric acid & then warm it $60^\circ - 70^\circ C$ and



This temp means one can hold the neck of conical flask by bare hand if temp is higher than this then one can't

hold the conical flask, the oxalic acid is heated to faster the reaction. At room temp. titration will be not be complete since reaction of $KMnO_4$ & oxalic acid will be very slow.

At $60^\circ - 70^\circ C$ temp conical flask is removed from the bunsen burner and initial reading of $KMnO_4$ from burette is already written in the copy hence immediately titration of oxalic acid with $KMnO_4$ is done. The titration is required to be completed in few mins so that solution remains hot. Initially the oxidation-reduction is balanced but as it is complete

no. of moles of O.A. is balanced by no. of moles of KMnO_4 , equivalence point is reached but we can't deduct with eyes. So the end point is performed by adding a slight excess of KMnO_4 that solution becomes pink which one can easily visualize with naked eye. A single excess drop of KMnO_4 is sufficient for observing the end point. Repeat the process to get three concurrent readings.

Titration of KMnO_4 vs std. oxalic Acid of strength 0.02476 N

S.No	IR. KMnO_4	FR.	Vol. of KMnO_4 used
1	0	11.0	11.0
2	0	11.0	11.0
3	0	11.0	11.0

$$(N_1 V_1)_{\text{Oxalic Acid}} = (N_2 V_2)_{\text{KMnO}_4}$$

$$0.02476 \times 10 = N_2 \times 11.0$$

$$0.02476 \times 10 = N_2$$

11

$$0.022 \text{ N} = N_2 \text{ (Normality of } \text{KMnO}_4 \text{)}$$

Step 3 Titration of Standardized KMnO_4 VS unknown oxalic Acid.

Take 10ml of unknown strength oxalic acid solution and add 1 Test Tube of H_2SO_4 . Heat the contents upto $60^\circ - 70^\circ\text{C}$ and then titrate with KMnO_4 present in the burette.

The end point is detected by observing pink color. Then the solution is thrown into the sink, wash the conical flask with water, and then take 10ml of unknown oxalic acid + H_2SO_4 and repeat with KMnO_4 . Repeat the procedure till three concurrent readings are obtained.

Sl. No.	IR.	FR	Vol. of standardized KMnO_4
1	0.0	12.0	12.0
2	0.0	12.0	12.0
3	0.0	12.0	12.0

$$(N_2 V_2')_{\text{KMnO}_4} = (N_3 V_3)_{\text{unknown oxalic acid}}$$

$$0.022 \times 12 = N_3 \times 10$$

$$\frac{0.022 \times 12}{10} = N_3$$

$$0.0022 \times 12 = N_3$$

$$0.0264 \text{ N} = N_3$$

Normality of ^{unknown} oxalic Acid solution = 0.0264
strength of unknown O.A. soln is

$$= 0.0264 \times 63 \text{ g/l}$$

$$= 1.6632 \text{ g/l}$$

Result:- The strength of unknown oxalic acid is 1.6632 g/l